

# PYCKO SCIENTIFIC LIMITED

31 London Road, Grantham, NG31 6EX

Telephone 01476 401992

bill@pycko.co.uk

www.pycko.co.uk

Your Alternative  
To The Obvious

## **MODEL 180-28**

### **SAMPLE HOLDER FOR MODEL 26-SERIES FRISKERS**



#### **Specifications**

INDICATED USE: repeatable geometry counting of wipes, filter paper, or slides

SAMPLE SIZE: 4.7 cm (1.8 in.) diameter maximum

SAMPLE TRAY: selectable spacing of 0.32 cm (0.125 in.), 0.64 cm (0.25 in.), 1.3 cm (0.5 in.), 2.5 cm (1 in.), and 5.1 cm (2 in.) from detector

SUGGESTED DETECTOR: Model 26 Integrated Frisker

CONSTRUCTION: anodized aluminium frame and sample tray

SIZE: 9.1 x 7.6 x 15.2 cm (3.6 x 3 x 6 in.) (H x W x L)

WEIGHT: 0.40 kg (0.9 lb)

#### **Model 26 Integrated Frisker**

The Model 180-28 Sample Holder is designed to accommodate the design of the Model 26 Integrated Frisker to make testing for contamination on wipes, filter paper, or slides simple and effective.

